



Lodging Program

Product Reference Guide

ECOLAB[®]

Everywhere It Matters.™



Delivering the kind of clean that's good for business.

Guest satisfaction and operational efficiency are vitally important to your success – and Ecolab can help you raise the bar for both. Our comprehensive cleaning and disinfecting solutions, combined with unparalleled service and expert training, makes Ecolab your go-to resource for all of your cleaning needs. Explore our wide-ranging solutions, specifically designed for the lodging business.



On-Premise Laundry Solutions 4

We help you maintain clean, soft linens and manage the total costs of your operation.



Housekeeping Solutions 15

We offer a complete program to meet all of your housekeeping needs.




Floor Care 22

We offer solutions that help you keep your floors gleaming all day.



Hand Hygiene 26

Our hand hygiene system combines innovative chemistry with effective training to help develop healthy habits.



Pool & Spa Solutions 31







Our solutions help you maintain clear, sparkling water and safe, stable chemistry.



Water Quality Management Solutions 33

Our filtration and softening solutions are designed for maximum safety and efficiency.

Aquanomic™ Solids

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Low Temp Laundry Solid Detergent	NPE & Phosphorus-free concentrated solid detergent is formulated to deliver premium results at a lower wash temperature, leaving linens whiter and brighter.	6100307	4-9 lb
 Low Temp Laundry Solid Destainer	Solid destainer formulated to retain bleaching power. For use with other Aquanomic low temperature products. Provides excellent bleaching and preserves linen strength and integrity. Helps lower overall linen replacement costs.	6100308	2-4 lb
 Low Temp Laundry Solid Clearly Soft™	Premium performance solid fabric softener with advanced chemistry. Delivers premium softness and long-lasting fresh fragrance. Integrated souring capability simplifies laundry operations and reduces the need for a separate sour in most applications.	6100309	2-6 lb
 Low Temp Laundry Solid Sour (3rd Product)	Solid sour for use with Aquanomic low temperature detergent and builder brings linen to a neutral pH. The phosphorus-free blend is effective for a wide range of conditions and on all fabrics.	6100311	2-6 lb
 Low Temp Laundry Solid Break	Solid alkaline break is effective at lower temperatures, tough on stains and gentle on linens. This concentrated formula inhibits scale formation and controls water hardness to ensure linens remain white over time.	6100310	4-6 lb
 Low Temp Laundry Solid Sour (4th Product)	Solid sour for use with Aquanomic low temperature detergent and builder brings linens to a neutral pH. This phosphorus-free blend is effective for a wide range of conditions and on all fabrics. For use in the optional fourth product dispenser if additional souring is required due to local water conditions.	6100312	4-6 lb

Aquanomic™ Liquids

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Low Temp Laundry Detergent	NPE & Phosphorus-free concentrated liquid detergent is formulated to deliver premium results at a lower wash temperature, leaving linens whiter and brighter.	6100445 6100446 (US only)	1-15 gal 1-55 gal
 Low Temp Laundry Destainer	Destainer formulated to retain bleaching power. For use with other Aquanomic low temperature products. Provides excellent bleaching and preserves linen strength and integrity. Helps lower overall linen replacement costs.	6100585 6100586 (US only)	1-15 gal 1-55 gal
 Low Temp Clearly Soft™	Premium performance liquid fabric softener with advanced chemistry. Delivers premium fabric softness and a fresh fragrance under tough water conditions.	6100588 6100589 (US only)	1-15 gal 1-55 gal
 Low Temp Laundry Sour	Liquid sour for use with low temperature detergent and builder brings linen to a neutral pH. The phosphorus-free blend is effective for a wide range of conditions and on all fabrics.	6100438 6100439 (US only)	1-15 gal 1-55 gal
 Low Temp Laundry Builder	Liquid alkaline builder is effective at lower temperatures, tough on stains and gentle on linens. This concentrated formula inhibits scale formation and controls water hardness to ensure linens remain white over time.	6100441 6100442 (US only)	1-15 gal 1-55 gal
 Low Temp Iron Controlling Sour	Liquid sour specially formulated for use in areas with high levels of iron in the water. Use with Aquanomic low temperature detergent and builder to bring linen to a neutral pH.	6100667 6100668 (US only)	1-15 gal 1-55 gal



Tri-Star™

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Tri-Star™ Detergent Oz	Concentrated, low odor, pH sensitive surfactant blend. Designed for use with ozone washing systems using low water temperatures.	16561	1-5 gal
 Tri-Star™ Duo	Liquid one-shot laundry detergent.	10395	1-2.5 gal
 Tri-Star™ Enzymatic Special	Enzymatic laundry detergent ideal for removing protein-based stains.	6100914	1-5 gal
 Tri-Star™ Flexylite™	Low alkaline detergent specially formulated for lightly soiled linens.	10356 16326	1-2.5 gal 1-5 gal
 Tri-Star™ L-2000 XP	High performance, alkaline, one-shot emulsion detergent with nonionic surfactants and conditioners. Performs in a wide range of wash temperatures and in challenging water conditions and is effective on both oily and particulate soils.	6100031	1-5 gal
 Laundry Neutral Detergent Plus	Neutral liquid laundry detergent penetrates and lifts off oily soils, whitens and brightens all fabrics and is effective in a wide range of temperatures.	6100172	1-5 gal
 Tri-Star™ Solar Brite™	High alkaline emulsion laundry detergent with nonionic surfactants, water conditioners and optical brighteners. Effective on all linen types and soils across a broad range of wash temperatures and water conditions.	6100106	1-5 gal
 Tri-Star™ Surlite™	A unique one-shot, concentrated laundry detergent with a patented blend of nonionic surfactants, the unique ratio of alkalinity to surfactant allow the suds and bleach steps to be combined on light and medium formulas, saving water and energy.	6100267 6100268 (US only)	1-2.5 gal 1-5 gal
 Tri-Star™ Break III Xtra	Concentrated alkaline builder with water conditioners. Inhibits scale formation and controls water hardness and yellowing from iron. (US only)	16558	1-5 gal
 Tri-Star™ Turbo Boost	A blend of natural solvent and surfactants without an aromatic smell. Provides an extra punch on problem, extra heavy oily or greasy soils.	14330	1-5 gal
 Tri-Star™ Laundri Destainer	Concentrated chlorinated bleach provides fast bleaching action, and removes difficult stains. Restores natural whiteness to fabrics and eliminates dull and dingy build-up. (US only)	15859 15982	4-1 gal 1-5 gal







Tri-Star™ (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Tri-Star™ Chlorlite	Stabilized hypochlorite-based stain remover for white linen. Will work in all water temperatures, but performance will be compromised below 120° F.	10261	1-2.5 gal
 Tri-Star™ Oxy-Brite™	Oxygen bleach efficiently removes stains on all white and colored washable fabrics. Helps prolong life by minimizing tensile strength loss due to chlorine retention. (US only)	16527	1-5 gal
 Tri-Star™ Neutralizer	Super concentrated laundry sour that adjusts linen pH to eliminate skin irritation. Offers improved iron control and removes yellowing for whiter whites without oversouring.	14893 12083	1-2.5 gal 1-5 gal
 Tri-Star™ Sour VII	Premier laundry souring agent for controlling high iron and hard water conditions. Use to help solve result issues from high iron or hard water conditions. Very aggressive acid.	16878	1-5 gal
 Tri-Star™ Advacare 120 Sanitizer/Sour	EPA-registered sanitizer and sour that provides high neutralizing capacity. (US only)	14299	1-4 gal
 Tri-Star™ Aqua Soft™	Cationic fabric softener and neutralizer, soluble in hot or cold water, contains optical brighteners and leaves a floral scent.	14895 12081	1-2.5 gal 1-5 gal
 Tri-Star™ Clearly Soft™	Liquid fabric softener that gives fabrics a luxurious softness and a fresh, clean springtime scent. Formulated to dose at higher levels, keeping linens their whitest with maximum absorbency.	14994 14995	1-2.5 gal 1-5 gal
 Tri-Star™ So Fresh™	Softener leaves linens soft and static free. Adds a gentle, clean and fresh fragrance.	14310 14320	1-5 gal 1-15 gal
 Laundri Bacstat Liquid Soft	Fabric softener concentrate, with residual self-sanitizer reduces odor causing bacteria under conditions of high relative humidity or wet contamination. Softens clothes, removes static cling and helps improve ironing time. EPA-registered. (US only)	17418	1-5 gal
 Polycrisp	General application laundry starch that imparts more body and a crisper feel to linens. Great on food and beverage linens at restaurants and hotels.	14411	1-5 gal

Tri-Star™ (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Natural Satin Starch	Natural laundry starch provides excellent finish to linen. Washes out more completely than synthetic sizings, preventing permanent buildup. (US only)	18838	1-5 gal
 Water Conditioner	Sequestering agent to control water hardness, iron and related heavy metals.	10400	1-5 gal
 Water Treatment MC	Laundry additive for high metal water conditions, effective in maintaining fabric whiteness and brightness, including removing sunscreen staining.	6100145	1-5 gal







Eco-Star™

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Eco-Star™ Detergent MP	30% active, multi-purpose detergent that contains strong wetting agents that penetrate the fabric to loosen soils for removal. The choice for most general applications where a less concentrated detergent is desired. (US only)	16124 16123 16122	1-15 gal 1-30 gal 1-55 gal
 Eco-Star™ Detergent Oz	Concentrated, low odor, pH sensitive surfactant blend. Designed for use with ozone washing systems using low water temperatures. (US only)	16175 16174	1-15 gal 1-55 gal
 Eco-Care™ Enzymatic Special	Concentrated enzyme detergent additive ideal for removing protein-based stains. (US only)	6100915	1-15 gal
 Eco-Star™ Flexylite™	Low alkaline detergent specially formulated for lightly soiled linens. (US only)	16325	1-15 gal
 Eco-Star™ L-2000 XP	High performance, alkaline, one-shot emulsion detergent with nonionic surfactants and conditioners. Performs in a wide range of wash temperatures, in challenging water conditions and is effective on both oily and particulate soils. (US only)	6100032 6100033	1-15 gal 1-55 gal
 Eco-Care™ LT Detergent	Concentrated detergent comprised of a special blend of surfactants, enzymes and oxygenated solvents. For optimal performance, use in lower temperature (120° F) wash conditions. (US only)	16195 16194	1-15 gal 1-55 gal




Eco-Star™ (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Laundry Neutral Detergent Plus	Neutral liquid laundry detergent penetrates and lifts off oily soils, whitens and brightens all fabrics and is effective in a wide range of temperatures. (US only)	6100173 6100171	1-15 gal 1-55 gal
 Eco-Star™ Solar Brite™	High alkaline emulsion laundry detergent with nonionic surfactants, water conditioners and optical brighteners. Effective on all linen types and soils across a broad range of wash temperatures and water conditions.	6100105	1-15 gal
 Eco-Star™ Builder C	Concentrated alkaline builder with advanced water conditioning polymer technology.	10430 10429	1-15 gal 1-55 gal
 Eco-Star™ Turbo Boost	A blend of natural solvent and surfactants without an aromatic smell. Provides an extra punch on problem, extra heavy oily or greasy soils.	16113 (US only) 16166	1-15 gal 1-55 gal
 Eco-Star™ Destainer	Concentrated chlorinated bleach provides fast bleaching action, and removes difficult stains. Restores natural whiteness to fabrics and eliminates dull and dingy build-up.	16146 16144	1-15 gal 1-55 gal
 Eco-Star™ Oxy-Brite™	Oxygen bleach efficiently removes stains on all white and colored washable fabrics. Helps prolong life by minimizing tensile strength loss due to chlorine retention. (US only)	16006 16148	1-15 gal 1-55 gal
 Eco-Star™ Neutralizer	Super concentrated laundry sour that adjusts linen pH to eliminate skin irritation. Offers improved iron control and removes yellowing for whiter whites without oversouring.	14889	1-15 gal
 Eco-Star™ Sour	High neutralizing capacity sour for the most economical use. Adjusts pH of garments to minimize potential of skin irritation and fabric browning. Should be used with care as overuse could burn linen. (US only)	16159 16158	1-15 gal 1-55 gal
 Eco-Star™ Sour VII	Premier laundry souring agent for controlling high iron and hard water conditions. Use to help solve result issues from high iron or hard water conditions. Very aggressive acid. (US only)	16659	1-15 gal
 Eco-Care™ Sour HC	Very aromatic mild sour. Adjusts the pH of garments to that of skin, helping to minimize potential skin irritation. Use with low alkaline detergents/water. A compatible, gentle sour for laundry equipment. (US only)	16179 16178	1-15 gal 1-55 gal

Eco-Star™ (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Eco-Star™ Advacare 120 Sanitizer/Sour	EPA-registered sanitizer and sour that provides high neutralizing capacity. (US only)	16104	1-12 gal
 Eco-Star™ Clearly Soft™	Liquid fabric softener that gives fabrics a luxurious softness and a fresh, clean springtime scent. Formulated to dose at higher levels, keeping linens their whitest with maximum absorbency.	14996 14997	1-15 gal 1-55 gal
 So Fresh™	Softener leaves linens soft and static free. Adds a gentle, clean and fresh fragrance.	14320	1-15 gal
 Eco-Star™ Polycrisp	General application laundry starch that imparts more body and a crisper feel to linens. Great on food and beverage linens at restaurants and hotels. (US only)	14412	1-15 gal
 Water Conditioner	Sequestering agent to control water hardness, iron and related heavy metals. (US only)	10401 10399	1-15 gal 1-55 gal
 Water Treatment MC	Laundry additive for high metal water conditions, effective in maintaining fabric whiteness and brightness, including removing sunscreen staining. (US only)	6100148	1-15 gal






Navigator™

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Solid Dual Surge	One-shot solid detergent has optimized surfactant/alkali ratio and contains a mixture of surfactants, which deliver high performance on most soil types, especially on oily soils.	16309	4-8 lb
 Solid Surge Plus	Ideal for hospitality or healthcare linens, this premium one-shot detergent combines alkalinity, water conditioners, a wetting agent, fabric brightening system and a light fragrance. Green Seal Certified.	17905	4-9 lb
 Solid Break Plus	Solid detergent booster provides an extra punch of higher alkalinity to wash wheel. (US only)	16009	4-6 lb




Navigator™ (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Solid Stain-A-Way™	Solid destainer with no spills. Solid formula maintains bleaching power longer than liquid destainers.	13649	2-4 lb
 Solid Oxy-Brite™	Solid color-safe oxygen laundry bleach. Green Seal Certified.	16016	2-3 lb
 Solid Navisour	This problem solving laundry sour contains a unique blend of powerful ingredients that provide maximum protection against yellowing of linen in high iron water conditions. Green Seal Certified.	16005	2-6 lb
 Solid Navisoft	Ideal for hospitality or healthcare linens, this highly concentrated mixture of fabric softener and sour leaves linens with a soft hand. Doesn't spot or stain like some liquid softeners.	16001	2-6 lb

Laundry Powders

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Laundri Brite	One-shot powder laundry detergent. Designed for front loading, agitation style and home-style machines. (US only)	10390	1-45 lb
 Laundri Perma-Brite Plus	Heavy duty, powder detergent and color safe bleach system for commercial and home-style machines. (US only)	10391	1-45 lb
 Royal Brite™ Plus	Powder enzyme laundry detergent with oxygen bleach. Triple enzyme system for powerful cleaning and stain removal. Color safe, faster to rinse and reduces sour costs.	6100917	1-45 lb
 Laundri Break	Powdered break additive boosts available alkalinity. Supercharged formula effective on grease, fat, protein soils and "problem" loads. Prevents discoloration and dulling films. (US only)	17459	1-45 lb
 Laundri Dry Bleach	Powdered organic bleach with fast dissolving crystals. Removes stubborn stains and helps eliminate yellow and graying.	17707	1-45 lb

Laundry Powders (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Laundri Oxygen Bleach	Powdered laundry oxygen bleach for colored linens.	10410	1-45 lb
 Laundri Soft	Concentrated powdered softener. Gives a soft, luxurious feel to fabrics while eliminating static cling, making ironing and folding operations faster. (US only)	17756	1-45 lb
 Laundri Starch	General purpose powder laundry starch for a fresh, crisp feel on natural fabrics or blends. Adds body to cotton uniforms and other laundry items, makes ironing easier.	18812	1-35 lb

Home Style Solid

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Brite White	Laundry detergent packaged in pre-measured, water-soluble packets to help eliminate waste and spilling. Enzyme stain removing formula and contains an odor counteractant to eliminate stale, musty odors.	61502222	250-1.2 oz
 Home Style Solid Laundry Detergent with Enzymes	High surfactant, solid laundry detergent for home-style washers. Effective against protein based soils.	6100150	4-4 lb
 Solid Destainer	Solid bleach destainer for homestyle washers.	14508	2-2.6 lb
 Solid Oxy-Guard	Offers the same safety and stability as Solid Destainer, plus the added protection of color safe destaining. Can be used to brighten colors at low level usage or for destaining colored linen at higher levels.	16021	2-2 lb







Stain Management

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 StainBlaster™ Starter Kit	Wall mounted kit contains one bottle of StainBlaster Enzyme Boost, StainBlaster Grease Remover, StainBlaster Makeup Remover and StainBlaster Destainer laundry pre-spotters.	6101110	1 kit
 StainBlaster™ Enzyme Boost	Laundry pre-spotter for common soils such as blood, grass, lipstick, food soils, cuffs & collars, etc.	6101068 6101069	4-22 oz 1-5 qt
 StainBlaster™ Grease Remover	Laundry pre-spotter for greasy soils such as tar, crayon, motor oil, ink, pizza, pasta sauce, etc.	6101070 6101071	4-22 oz 1-5 qt
 StainBlaster™ Makeup Remover	Laundry pre-spotter for foundational makeup, sunscreen, medicine, baseball clay, etc.	6101072 6101073	4-22 oz 1-5 qt
 StainBlaster™ Destainer	Laundry pre-spotter for everyday bleachable stains such as coffee, tea, mildew, wine, etc.	6100374 6100395	4-22 oz 1-5 qt
 StainBlaster™ Rust Remover	Laundry pre-spotter for use on rust or iron stains. (Item not included in starter kit)	6100398	4-22 oz
 StainBlaster™ Multi-Purpose	Laundry pre-spotter for use on food grease, oil, body soils, lipstick, etc. (Item not included in starter kit)	6100370 6100371	4-22 oz 1-5 qt
 StainBlaster™ Power Pak Plus Reclaim White	Laundry reclaim solution for white linen. Pre-measured "Ready to Use" packets for proper dosing.	6100909	12-1.2 lb
 StainBlaster™ Power Pak Plus Reclaim Color	Laundry reclaim solution for colored linen. Pre-measured "Ready to Use" packets for proper dosing.	6100910	12-1 lb
 StainBlaster™ Power Pak Reclaim Rust Remover	Laundry reclaim solution for linen with rust or iron stains. Pre-measured "Ready to Use" packets for proper dosing.	6100911	12-0.5 lb



Dryer

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 X-Static™ Dryer Blocks	Powerful static control dryer blocks reduce static, leave a fresh, clean fragrance in the linen and enhance softness.		
	Mountain Fresh (12-pack blocks)	10311	12-150 g
	Mountain Fresh (6-pack blocks)	10309	6-150 g
	Neutral Scent (12-pack blocks)	10313	12-150 g
	Neutral Scent (6-pack blocks)	10312	6-150 g
	Woody Floral (12-pack blocks)	10317	12-150 g
	Woody Floral (6-pack blocks)	10316	6-150 g
	X-Static Dryer Block Dispenser	92592193	1 each




Laundry Dispensers

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Aquasonic Solids Dispenser	Aquasonic branded dispenser	9259-1101	1 unit
 Tri-Star™ Liquids Dispenser	ALCS 4 pump, 110V	9258-1007	1 unit
 Eco-Star™ Liquids Dispenser	Eco-Star and Aquasonic Liquids	9259-2427	1 unit
 Navigator™ Dispenser	Dispenses 3 solid laundry products	9259-1020	1 unit
 Home-Style Solid Dispensers – Solid System 1	Dispenses 1 home-style solid product	9253-1193	1 unit
 Home-Style Solid Dispensers – Solid System 2	Dispenses 2 home-style solid products	9253-1177	1 unit




Bathroom Disinfectants

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 66 Heavy Duty Alkaline Bathroom Cleaner and Disinfectant	Heavy duty alkaline cleaner with a pleasant floral fragrance for bathroom surfaces. Germicidal detergent removes common soils and soap scum. Disinfects at 8-12 oz. per gallon. EPA-registered.	14914 10158 20587 (Canada only)	1-2.5 gal 2-2 L 2-2 L
 73 Disinfecting Acid Bathroom Cleaner	EPA-registered heavy-duty and hospital grade bathroom disinfectant ideal to remove hard water build-up and soap scum from showers. Disinfects in 5-minutes compared to industry standard of 10-minutes. Offers a positive user experience with a soapy fragrance. Ideal for cleaning tough bathroom soils at low concentrations for economical results.	6101041 6101040	1-2.5 gal 2-2 L

Bathroom Cleaners

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Acid Bathroom Cleaner	Formulated from plant-derived surfactants to provide exceptional cleaning throughout bathrooms. Effectively removes dirt, soap scum and hard water deposits. Safe for use on bathroom chrome, porcelain, vinyl, tile and fiberglass. Green Seal Certified.	6100292 6100293	1-2.5 gal 2-2 L
 Neutral Bathroom Cleaner	Bio-based daily bathroom cleaner removes soil from toilet, shower walls and floors. Provides excellent cleaning through a plant-derived surfactant. Green Seal Certified.	6100275 6100276	1-2.5 gal 2-2 L
 Oasis Pro™ 70 Marble-Safe Cleaner	Gentle cleaner, safe for sensitive stone surfaces such as marble, terrazzo and limestone. Easily removes soils typically found in bathroom areas. Pleasing clean fragrance.	21032	2-2 L

All Purpose Disinfectants

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Peroxide Multi-Surface Cleaner & Disinfectant	True multi-purpose solution designed to drive operational efficiencies. EPA-registered product has a comprehensive claim set with two dilution rates for disinfection: 4 oz. per gallon (for 5 minute kill time) and 6 oz. per gallon (for 3 minute kill time). Recommended for cleaning and disinfecting hard surfaces, including streak-free performance on glass.	6100791	2-2 L
		6100793	1-2 gal
		6100897 (Canada only)	2-2 L
		6100894 (Canada only)	1-2 gal
 14 Antibacterial All-Purpose Cleaner	Help reduce the spread of microorganisms on hard surfaces. Concentrated, one-step, disinfectant cleaner is effective against a broad spectrum of bacteria, is virucidal, fungicidal and inhibits the growth of mold and mildew and their odors when used as directed.	14653	1-2.5 gal
		14205	2-2 L
		6100116 (Canada only)	2-2 L
 Neutral Disinfectant Cleaner	Multi-purpose, neutral pH, germicidal detergent and deodorant effective in hard water. Disinfects, cleans and deodorizes in one labor-saving step. Effective on hard surfaces and can also be used on non-porous athletic mats, exercise equipment and personal safety equipment (helmets, etc). Leaves a clean fresh scent. EPA-registered.	14542	1-2.5 gal
		6100388 (Canada only)	1-2.5 gal





All Purpose & Glass Cleaners

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Peroxide Multi-Surface Cleaner	Daily multi-purpose cleaner concentrate offers excellent cleaning and leaves a streak-free shine through the use of peroxide and a plant-based surfactant. Formulated to clean windows, glass, countertops, glazed porcelain, metal fixtures, walls, impervious floors and other hard surface materials. Product is not registered to sanitize or disinfect.	6100302	1-2.5 gal
		6100925	2-2 L
 Orange Force™	Streak-free, multi-surface cleaner with a pleasant orange fragrance. Removes tough greasy soils and films from multiple surfaces.	14559 (137)	1-2.5 gal
		21038 (16)	2-2 L
 Glass Cleaner	Glass cleaner concentrate formulated from plant-derived surfactants effectively removes tough soils including smoke and grease films. Streak-free, fast-drying performance on windows, mirrors, glass, plastic, acrylics and painted surfaces. USDA Certified BioBased product. Green Seal Certified.	6100288	1-2.5 gal
		6100289	2-2 L



Odor Control

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE		
 ZephAir Air & Fabric Refreshers	ZephAir Air and Fabric Refreshers tackle the toughest malodors, especially smoke and food odors, eliminating them quickly, safely and completely. Won't leave a sticky residue on fabrics and surfaces. <i>Clean White Cotton: The tool for tackling food and fitness room odors.</i> <i>Mountain Mist: Specializes in eliminating smoke and pet odors with a fresh, brisk and clean fragrance.</i>				
		Oasis Pro™ ZephAir Clean White Cotton	6100133	2-2 L	
		Oasis Pro™ 55 ZephAir Mountain Mist	6100984	2-2 L	
		Oasis™ ZephAir Clean White Cotton	6100132	1-2.5 gal	
		Oasis™ 271 ZephAir Mountain Mist	6100985	1-2.5 gal	
		ZephAir Clean White Cotton Ready-to-Use	6100129	6-32 oz	
 Room Refreshers	Specifically designed to neutralize odors while providing a pleasant and subtle, yet lasting clean scent.				
		ZephAir Mountain Mist Ready-to-Use	12043	6-32 oz	
		Oasis™ Clean Escape	6100830 (US only)	1-2.5 gal	
		Oasis™ Laundry Fresh Room Refresher	6100142	1-2.5 gal	
		Oasis™ Lavender Paradise Room Refresher	6100143 (US only)	1-2.5 gal	
		Oasis™ Morning Breeze Room Refresher	6100141	1-2.5 gal	
		Oasis™ Relaxing Spa	6100828 (US only)	1-2.5 gal	
		Oasis Pro™ Clean Escape Room Refresher	6100831	2-2 L	
		Oasis Pro™ Laundry Fresh Room Refresher	6100136	2-2 L	
		Oasis Pro™ Lavender Paradise Room Refresher	6100137	2-2 L	
		Oasis Pro™ Morning Breeze Room Refresher	6100135	2-2 L	
		Oasis Pro™ Relaxing Spa	6100829	2-2 L	
		 First Impression™ Micro Aerosol Odor Neutralizer	Aerosol technology neutralizes malodors in hallways and public restrooms. Programmable, automated dispensing provides continuous odor control with a wide variety of fragrances.		
				Micro Aerosol Mandarin Burst w/Dispenser	30074
Micro Aerosol Mandarin Burst	30078			12-1.8 oz	
Micro Aerosol Spring Showers w/Dispenser	30073			1 kit	
Micro Aerosol Spring Showers	30077			12-1.8 oz	
Micro Aerosol Summer Linen w/Dispenser	30071			1 kit	
Micro Aerosol Summer Linen	30076			12-1.8 oz	
Micro Aerosol Variety Pack	15093			12-1.8 oz	



Ready-to-Use Disinfectants

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Asepti-Wipe II	Broad-spectrum pre-moistened wipes. A convenient method for disinfecting hard surfaces. Good for use on remote controls and telephones. EPA-registered.	63306725 61027330	12-180 ct 6-160 ct
 TB Disinfectant Cleaner Ready-to-Use	Virucidal disinfectant with a broad set of claims. Formulated to disinfect and deodorize hard, non-porous, environmental surfaces. Effective against Norwalk virus (30 seconds), HBV, TB and many other organisms as listed on the label. EPA-registered. (US only)	65443556	12-32 oz
 CBC™ Plus Bowl Cleaner	Mild acid bowl and disinfecting cleaner that combines detergent and wetting agents coupled with a phosphoric acid base. Eliminates rust, lime, organic deposits, and offensive odors from toilet bowls, urinals, wash basins, bathtubs, shower stalls and drinking fountains. Contains corrosion inhibitors.	61143699 5705 (Canada only)	12-32 oz 12-32 oz
 Thickened Non-Acid Bowl/Bathroom Cleaner	Powerful, thick formula that clings to vertical surfaces for improved cleaning and disinfecting. Safe for use on synthetic surfaces, chrome fixtures and for general cleaning and disinfecting of bath and shower areas including toilet bowl. Will not damage carpeting.	61143665 5704 (Canada only)	12-32 oz 12-32 oz

Ready-to-Use Bathroom Cleaners

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Aim Encounter Urinal Block with Screen	Releases cleaning agents with each flush, inhibits deposit buildups and uses odor counteractant technology to eliminate odors. Formulated for longer life.	61144200	12-100g
 Gel Pro™	Powerful phosphorus-free gel bathroom cleaner with a fresh lemon fragrance and mild mildew stain remover has more destaining power than many consumer products. Clinging gel formula provides excellent surface contact and effectively cleans vertical tile and grout surfaces, no-slip strips and mats, tubs, showers and caulking.	6100574	12-16 oz
 Lemon Lift	Heavy duty cleaner and destainer specially formulated with chlorine bleach to remove even the toughest mold and mildew stains. Lifts away soils and leaves a fresh, lemon scent. (Not a disinfectant)	11490	12-20 oz




Ready-to-Use Bathroom Cleaners (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Lemon-Eze™	Crème cleanser quickly removes difficult spots, soap scum build-up, rust, grease and more. Leaves surfaces clean and shiny and effectively cleans without harsh abrasives. Will not scratch or dull surfaces when properly used. Light, lemon fragrance.	13094	12-32 oz
 Bingo™ Liquid Drain Cleaner	Heavier-than-water formulation reaches and dissolves hair, protein, grease, oil and soap. Harmless to pipes and septic tanks. Use to open clogged pipes and perform routine maintenance on plumbing systems.	61035112	12-32 oz





Ready-to-Use All Purpose Cleaners

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Orange Force™ Ready-to-Use	Streak-free, multi-surface ready-to-use cleaner with a pleasant orange fragrance. Removes tough greasy soils and films from multiple surfaces. (US only)	15639	6-32 oz
 Power Force Heavy Duty Degreaser Cleaner	Heavy duty, premium degreaser that eliminates even the toughest soils, from crayon to ink stains. Uses a combination of solvents and surfactants to deliver premium cleaning performance. (US only)	14544	6-32 oz
 Behold Lemon Furniture Polish	Professional formula penetrates dust, oil and wax to remove stains and smudges. Protects surfaces and adds new luster to furniture.	61196406	6-16 oz
 Radiance Furniture Polish	High performance furniture polish provides deep cleaning and leaves a long lasting finish that repels dust. Special formula fills in nicks and micro-scratches and leaves behind a clean, fresh scent.	10120	12-17 oz
 Endust Dust Control	Makes a dust magnet of your cloth or mop, giving it three times its normal dust-holding capacity. No-wax formula is ideal for regular, cost-effective cleaning and polishing. Prevents dust scattering and resettling.	61196291	6-16 oz
 Java Klenz	Specially formulated, concentrated granular compound removes stubborn stains created by coffee and tea for a positive guest impression. Removes rancid coffee oils and tannins which can destroy the "fresh brewed" aroma and cause off-flavors. Safe on plastic, china and metal surfaces and effective in both hard and soft water conditions.	17585	100-0.25 oz

Housekeeping Dispensers

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Select 4 Dispensing System	Each dispenses up to four 2.5 gallon products to either spray bottle or mop bucket. Air Gap available.	92151149 92151150	1 unit 1 unit
 Ultra 1 Dispensing System	Each dispenses one 2.5 gallon product to either spray bottle or mop bucket. Air Gap available.	92151151 92151152	1 unit 1 unit
 Oasis Pro™ Integrated Spray Bottle Dispenser	Dispenses up to four 2 liter bags to spray bottles.	Oasis Pro™ Integrated AVB Dispenser 92151143 Oasis Pro™ Integrated Flex Gap Dispenser 92151142	1 unit 1 unit







Housekeeping Tools

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Housekeeping Tools Complete Tool Kit	Kit Contains: Extendable Pole, Swivel Pad, Foam Squeegee, Tape Roller, Non-Abrasive Scrubbing Pad, Microfiber cleaning pad, Power Eraser™ Plus pad, Pump-Up Foamer/Sprayer, 2 blue, 2 pink and 2 yellow Microfiber cloths, Laundry mesh bag, Ecolab Rinse Cup.	9621-1104	Kit
Housekeeping Tools Starter Kit	Tools only – Extendable Pole, Foam Squeegee, Adhesive Tape Roller, and Swivel Pad holder.	9621-1110	Kit
 Telescopic Pole	Easily extends from 2'-4' to help room attendants reach high walls and other difficult-to-clean areas. Attaches to tool heads with a specially designed quick connect fitment.	9621-2000	1 each
 Swivel Cleaning Pad	Tool head with 3 easy to change pads is ideal for cleaning tubs, walls and other bathroom surfaces. Comes with Non-Abrasive Scrubbing pad, Power Eraser Plus pad and Microfiber pad.	9621-1109	3 pack
 Foam Squeegee	14-inch foam squeegee is ideal for quickly drying bathtub surfaces, including tile. Quick connect fitment easily connects to the Extendable Pole, sold separately.	9621-1106	3 pack



Housekeeping Tools (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Adhesive Tape Roller	Use as a finishing touch for removal of any hairs or small particles, especially in the bathroom. Includes 3 rollers and 3 tape rolls. Quick connect fitment easily connects to the Extendable Pole, sold separately.	9621-1107	3 pack
 Pump-Up Foamer/Sprayer	A 32-oz pump-up foamer and sprayer in one unit. Reduces repetitive motion injury by reducing the need to pull a standard trigger sprayer.	9265-2602	1 each
 Microfiber Cloths	Microfiber weave improves performance of cleaners, and leaves a lint-free, streak-free surface. Comes in blue, yellow and pink to match product color coding.	Blue 9621-2290 Pink 9621-2292 Yellow 9621-2293	12 pack 12 pack 12 pack
 Power Eraser™	Mildly abrasive pad quickly and efficiently eliminates tough stains, marks and soap scum on hard surfaces, such as walls, doors, countertops, shower doors, tubs and fixtures. Also great for luggage marks on baseboards.	9290-2035 9297-2030	25 pk 12 pk
 Non-Abrasive Scrubbing Pad Refills	Ideal for cleaning tubs and shower walls. Use with quick connect swivel pad.	9621-1000	12 pack
 Power Eraser™ Plus Refills	Designed to scrub tough-to-clean stains on tub bottoms. Use with quick connect swivel pad.	9297-2030	12 pack
 Microfiber Cleaning Pad Refills	Great for cleaning bathroom floors. Use with quick connect swivel pad.	9297-2031	12 pack
 Adhesive Tape Refills	High quality sticky paper rolls with perforations to allow easy removal of dirty papers.	8999-0011	12 pack
 Ecolab Rinse Cup	16 oz. durable plastic cup made for rinsing surfaces, especially bathtubs, after cleaning.	9263-2026	1 each
 Laundry Mesh Bag	Mesh bag with drawstring closure. Ideal for laundering dirty microfiber cloths.	9298-3001	6 pack

Carpet Spotters

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Revitalize Miracle Spotter Starter Kit	Special twin cartridge design blends active components on the spot, providing results not available with single component chemistry. Ideal for the "unknown" spot. Use on all synthetic carpet types.	14511 14512	2-20 oz 4-20 (refills)
 Revitalize Spotting Kit	Easy-to-use, lightweight spotting kit. Hang on cleaning carts for easy transportation. Kit includes: Revitalize Multi-Purpose Spotter, Beverage and Healthcare Spotter, Gum Remover, Multi-Purpose Solvent Spotter and spotting guide.	14618	1 kit
 Revitalize Gum Remover	Effectively freezes gum and candle wax to be quickly chipped away. Leaves no residue. Non-combustible and non-butyl. Contains no CFCs. Use on all carpet types.	61495088	12-6 oz
 Revitalize Multi-Purpose Solvent Spotter	Powerful solvent spotter for heavy grease and oil deposits. Rinses easily with water and leaves no oily residue – no resoiling. Use on all synthetic carpet types.	10386	6-24 oz
 Revitalize Enzyme Spotter and Odor Eliminator	Highly-effective enzyme-based product for removal of organic stains and elimination of related odors. Great for elimination malodors at the source, such as urine, with no fragrance cover-up. Use on all synthetic carpet types.	61495092	6-24 oz
 Revitalize Rust Spotter	Effectively removes iron stains from carpet. Does not contain hydrofluoric acid. Use on all synthetic carpet types. (US only)	61495019	6-24 oz




Carpet Cleaners

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Revitalize Encapsulation Carpet Cleaner	Low-moisture interim cleaner that uses a polymer technology to encapsulate soil and prevent it from reattaching to the carpet. Soil is removed during the next scheduled vacuuming. Great for quick turnaround in high-traffic areas with a dry time of less than one hour. Use on all synthetic carpet types. For use with pump-up sprayer, encapsulate equipment and/or carpet brush.	61495099	4-1 gal
 Revitalize Carpet and Upholstery Prespray Cleaner	Effective spin bonnet cleaner for use with pump-up sprayer or conventional spin bonnet machines. Use as a pre-treatment spray on traffic lanes prior to extraction or shampooing. Use on all carpet types and upholstery. Use with pump-up sprayer followed by extraction. Green Seal Certified.	61495060 61495062 61495069	4-1 gal 2-2 L 1-2.5 gal





Carpet Cleaners (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Revitalize Carpet and Upholstery Extraction Cleaner	Superior extraction chemistry with excellent penetration, soil lifting and cleaning performance. Use on all carpet types and upholstery. Use with carpet extractors.	61495074 61495077 61495078	4-1 gal 2-2 L 1-2.5 gal
 Revitalize Carpet Shampoo	Rich-foaming carpet shampoo lifts embedded soils to the surface. Chemistry includes unique odor counteractants to freshen carpet. Use on all synthetic carpet types. Use with rotary or mistfoam machines.	61495122	4-1 gal
 Revitalize Perfusion Carpet and Upholstery Cleaner	Dual-action chemistry eliminates the need for pre-spotting when used as an extraction cleaner. Also effective when used as a pre-spray or spotter. It's portable, so extractors may be filled onsite without a dilution station. Not wool or wool-blend safe. For use with carpet extractors, pump-up sprayers and spray bottles.	10379	1-10 lb

Carpet Specialty Solutions

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Revitalize Anti-Soil Treatment	Specialized treatment for protecting carpets. Repels oil and water-soluble stains. Ready-to-use formula contains DuPont Teflon® and is non-combustible. Use on all synthetic carpet types. For use with pump-up sprayers.	61495025	4-1 gal
 Revitalize Defoamer	Instant defoaming action improves operational efficiency by reducing need to empty recovery tank. Use in carpet extractor recovery tanks.	61495113	4-1 gal
 Revitalize Non-Browning Treatment	Neutralizer rinse eliminates browning caused by water staining and over-wetting. Also effective to neutralize ice melt compounds and tracked in salts on carpets. Use on all synthetic carpet types. For use with carpet extractors or pump-up sprayers.	61495036	4-1 gal
 Revitalize Rug and Room Deodorizer	Dry, pleasant scented deodorizer neutralizes malodors caused by smoke, mustiness and general occupancy. Use on all carpet types. For use as a treatment prior to vacuuming.	61495141	12-22 oz







Floor Cleaners & Restorers

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 High Performance Neutral Floor Cleaner	Ultra low-foaming, neutral floor cleaner designed for daily use in automatic scrubbers and mop buckets. The powerful blend of surfactants can be used with hard or soft water on all floor surfaces and is designed to remove soils without leaving a film or haze. Cleans without removing floor finish or dulling gloss and does not require rinsing. Green Seal Certified.	6100035 6100036	2-2 L 1-2.5 gal
 Maxx Dual Action Floor Cleaner	Great for scrub and recoat applications. Ultra low-foaming chemistry increases auto scrubber efficiency.	6100428 6100659	1-2.5 gal 2-2 L
 Xpressway™ Neutral Floor Cleaner	Ready-to-use neutral floor cleaner for daily use on all floor types. Use with the Xpressway™ Mop. (9202-1080)	11222	6-1 L
 Prep-Pak™ Neutralizer/Conditioner	Promotes proper adhesion and performance of floor finish by removing alkali stripper residue before recoating. Also dissolves ice melt compound from floor surfaces and lime scale from glass, tile and metal.	61191402	180-0.5 oz
 Sanitizing Wash 'n Walk	EPA-registered cleaner and sanitizer provides total management of kitchen floors and drains in a single product. Enzyme based floor cleaner provides immediate cleaning of all soil types and extended cleaning of organic grease, fats and oils, thereby reducing risks of slips and falls. Sanitizer component of the formula kills 99.9% of pathogens, reducing cross contamination risk, malodor, etc.	6100731 6100734	2.5 gal 15 gal



Stone Cleaners

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Pinnacle™ Hard Surface Polish	Ready-to-use product to clean and polish all natural stone, glazed ceramic, laminates, plastics, vinyl, fiberglass, stainless steel and most painted surfaces. Provides a long-lasting shine, protects the surface and has a fresh citrus scent. Not to be used on floors.	62700138	12-18 oz
 Pinnacle™ Floor Cleaner/Conditioner	Helps develop a natural patina and increased gloss while providing added protection against soiling. Special conditioners enhance the natural color of marble, terrazzo and granite. For daily use. (US only)	62770217	4-1 gal
 Pinnacle™ Intensive Tile Cleaner	Alkaline formula removes stubborn soils and grease from ceramic tiles and grout. For interim maintenance.	62770211	4-1 gal

Stone Polishers

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Pinnacle™ Marble Polish Paste	Easy-to-use compound that restores polish of worn marble, travertine and other carbonate-based stone. Polishes, conditions and offers initial water and stain resistance without damaging stone surfaces.	62770418	1-10 lb
 Pinnacle™ Marble Polish Prep	Non-acid formula removes light scratches in marble when final polishing compounds are no longer sufficient. Requires follow-up with Marble Paste.	62770420	2-5 lb
 Pinnacle™ Marble & Terrazzo Crystallizer	Ideal for when polished marble and terrazzo begins to dull. Produces a high sheen and restores color while increasing the durability of the polish. Use as needed.	62771300	4-1 qt
 Pinnacle™ Granite Polish Cream	Easy-to-use, ready-mixed formula repairs the light scratches and wear patterns in traffic lanes. Restores lightly-worn granite while enhancing water and stain resistance.	62771433	4-1 qt
 Pinnacle™ Tile Polish #1	Step 1 of a 2-step process. Cleans and prepares tile for polishing. Also cleans the grout so the entire area being treated is prepped for Tile Polish 2. Use as needed.	62770411	1-10 lb
 Pinnacle™ Tile Polish #2	In conjunction with Tile Polish 1, this ready-to-use product adds a shiny, natural gloss to most unglazed ceramic, porcelain and quarry tile floors. Use as needed.	62770413	1-10 lb

Stone Specialty Solutions

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Pinnacle™ Grout Impregnator	Ideal for repelling moisture and protecting grout lines from oil, grease and soil. Impregnators absorb into the pores. Use as needed with Pinnacle™ Intensive Tile Cleaner.	62700162	4-1 gal
 Pinnacle™ Traction	Helps increase the floor's coefficient of friction. Use as needed.	62700104	4-1 gal


Nexa General Use Hand Soaps

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Clean & Smooth	Green Seal Certified lotion hand soap that provides rich, moisturizing lather to remove stubborn dirt while leaving hands fresh and soft.	6100706	6-750 mL
		6100707	4-1250 mL
 Foam Hand Soap	Green Seal Certified gentle foaming hand soap that is ideal for public restrooms.	6100708	6-750 mL
		6100709	4-1250 mL

Nexa Antibacterial Hand Soaps

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Antibacterial Clean & Smooth	High-quality antibacterial soap that is gentle on hands and effective against germs.	6100711	6-750 mL
		6100712	4-1250 mL
 Antibacterial Lotion Hand Soap	New and improved formula featuring chloroxylenol (PCMX) to deliver effective antimicrobial efficacy. It lathers quickly and rinses freely, minimizing residual soap left on the hands. (US only)	6100713	6-750 mL
		6100714	4-1250 mL
 Antibacterial Foam Hand Soap	Pre-lathered antibacterial soap with a light citrus scent that kills microorganisms.	6100715 (US only)	6-750 mL
		6100716 (US only)	4-1250 mL
		6100754 (Canada only)	6-750 mL
		6100755 (Canada only)	4-1250 mL
 Medicated Foam Hand Soap	Antibacterial hand wash that meets FDA requirements for Healthcare Personnel Handwash. Contains aloe vera and emollients to help moisten and soothe the skin. (US only)	61100717	6-750 mL
		6100718	4-1250 mL
 Gentle Antibacterial Foam Hand Soap	Protects against microorganisms in foodservice settings and features a citrus scent. (US only)	6100719	6-750 mL
		6100720	4-1250 mL


Nexa Hand Lotions

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Revitalizing Skin Lotion	Advanced formula designed to be non-greasy and quick-absorbing to provide superior moisturization. Fragrance-free.	6000062	6-750 mL



Nexa Hand Sanitizers

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Foam Hand Sanitizer	Formulated with 62% ethyl alcohol, promotes skin health, reduces the spread of infections and is a USDA-certified biobased product. (US only)	6100727	6-750 mL
		6100728	4-1200 mL
 Food Service Hand Sanitizer	FDA food code-compliant (non-alcohol based). Contains skin conditioners and emollients to help moisturize and improve skin health. Ideal for use in the kitchen or foodservice settings.	6100729 (US only)	6-750 mL
		6100730 (US only)	4-1250 mL
		6100758 (Canada only)	6-750 mL
		6100759 (Canada only)	4-1250 mL
 Instant Hand Sanitizer Gel	This formula is more than 99.5% effective in in-vitro testing against a broad range of foodborne illness-causing bacteria, helping to minimize the spread of contamination. It contains a unique ingredient for superior after-feel. (US only)	6100722	6-750 mL
		6100723	4-1200 mL
 Gentle Hand Sanitizer Gel	Formulated with 70% w/w ethyl alcohol, kills 99.99% of pathogens tested within 15 seconds. Built with vitamin E, moisturizers and five emollients to help maintain skin moisture. (US only)	6100725	6-750 mL
		6100726	4-1200 mL



Nexa Body Shampoos

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Total Body Shampoo	A protein-enriched, conditioning hair and body wash formulated to soothe and moisturize hair and skin. Can be diluted for use as a no-rinse bed/bath or perineal cleaner.	6000064	6-750 mL

Nexa Dispensers

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Touch-Free Dispenser	Energy saving dispenser requires fewer battery changeouts. Patented pump reduces clogs and drips and helps to avoid costly product waste.		
	Classic Touch-Free Dispenser (black)	9202-1192	1-1250 mL
	Classic Touch-Free Dispenser (white)	9202-1193	1-1250 mL
	Compact Touch-Free Dispenser (black)	9202-1190	1-750 mL
	Compact Touch-Free Dispenser (white)	9202-1191	1-750 mL
 Manual Dispenser	Patented pump pulls excess product back in, reducing clogs and drips. Antimicrobial protection on the push bar continuously inhibits growth of bacteria, mold & mildew.		
	Classic Manual Dispenser (black)	9202-3091	1-1250 mL
	Classic Manual Dispenser (white)	9202-3093	1-1250 mL
	Compact Manual Dispenser (black)	9202-3087	1-750 mL
	Compact Manual Dispenser (white)	9202-3088	1-750 mL

Nexa Dispenser Accessories

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
Drip Tray	Drip Tray (black)	9202-2989	1 unit
	Drip Tray (white)	9202-2993	1 unit
Lock and Key Kit	Lock and key kit.	9202-3113	2 keys, 10 locks
 Product Badges	Antibacterial Hand Soap	9202-3107	6 badges
	Hand Sanitizer	9202-3106	6 badges
 Sanitizer Station	Includes Nexa Classic Touch-Free Dispenser (white).	9202-6037	1 unit

Nexa Concentrates

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Nexa Concentrates System – 2L Bag	Concentrates provide sustainability and cost savings. System fills uniquely designed refillable bottles for use in Nexa dispensers.		
	Cosmetic and Sanitizer Antibacterial (US only)	9202-1200 9202-1201	1 unit 1 unit
 Refillable Bottles	Refillable bottles feature a sealed system and in-chemistry preservatives to ensure safe and effective application time-after-time.	9202-3330	4-1250 mL
		9202-3329	10-1250 mL
		9202-3328	4-750 mL
		9202-3327	10-750 mL
 Concentrated Foam Hand Soap	Concentrated general use foam soap that can be diluted up to a 10:1 dilution rate.	6100740	2-2 L
	 Concentrated Antibacterial Foam Hand Soap	Concentrated FDA-regulated foaming antibacterial soap with a benzalkonium chloride active and a light citrus scent. It can be diluted up to a 5:1 dilution rate. (US only)	6100867
 Concentrated Foam Hand Sanitizer		FDA food code-compliant and alcohol-free concentrated foam hand sanitizer with a benzalkonium chloride active. It can be diluted up to a 10:1 dilution rate. (US only)	6100868

Counter Mount – Dispensers & Soaps

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Counter Mount Dispenser	Touch-free dispenser for counter mount platform.	9263-1012	1 unit
 Antibacterial Foam Hand Soap	Antibacterial foam hand soap moisturizes and reduces bacteria. Contains moisturizers.	6100202	2-1600 mL
		6100272 (Canada only)	2-1600 mL
 Foam Hand Soap	Green Seal Certified gentle-foaming, cosmetic hand soap.	6100122	2-1600 mL

Manual Bulk – Dispensers & Soaps


PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Bulk Dispenser	Dispenser for Manual Bulk Refillable Liquid Hand Soaps.	9263-1970	1 unit
 Clean and Smooth Liquid Hand Soap	Gentle and creamy lotion soap enriched with emollients and moisturizers. Green Seal Certified.	15875	4-1 gal
 Antibacterial Clean and Smooth	Antibacterial liquid hand soap moisturizes and reduces bacteria.	14886	4-1 gal

Additional Sanitizers

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Digisan™ Aerosol Foam Hand Sanitizer	Quick foam hand sanitizer can be wall mounted or stand alone; great for tight spaces. (US only)	24832	12-7 oz
 Gel Hand Sanitizers	Non-dispensed waterless hand sanitizer with moisturizers to condition and protect. (US only)	6030352 6000003 6000004	144-1.2 oz 72-4 oz 12-540 mL

Aqua Balance™ Products

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Aqua Balance™ Total Alkalinity	Bring pool water alkalinity to the desired range of 80-120 ppm, helping to maintain pH and chlorine levels and keep water balanced, clean and safe for guests.	17415 17416 17891 (Canada only)	3-10 lb 1-25 lb 3-10 lb
 Aqua Balance™ pH Plus	Bring pool water pH to the ideal range (7.2-7.8) to keep swimmers comfortable and equipment lasting longer.	17401 17402 17421 (Canada only)	1-25 lb 3-10 lb 3-10 lb
 Aqua Balance™ Cal Hypo Shock	Superchlorinate (or "breakpoint chlorinate") pool water when combined chlorine levels are too high. (EPA Registration #7124-29-1677)	17413 14623 (Canada only)	30-1 lb 30-1 lb
 Aqua Balance™ Dechlor Reducing Agent	Quickly reduces the chlorine level in a 10,000 gallon pool by 16 ppm with one pound of product.	17411 17412 17422 (Canada only)	6-3 lb 6-10 lb 3-10 lb
 Aqua Balance™ Defoamer	Effectively break down oils to reduce foam and improve spa's appearance.	17403	4-16 oz
 Aqua Balance™ Basic Hardness	Protect pool from etching when soft water readings fall below the ideal 200-400 ppm range of calcium hardness. (US only)	17409 17410	3-10 lb 1-25 lb
 Aqua Balance™ Chlorine Free Oxidizer	Removes contaminants and helps eliminate chlorine odor. Swim shortly after use.	17414 17419 (Canada only)	12-1 lb 12-1 lb
 Aqua Balance™ Conquest Stain Remover	Removes stains and clears mineral-discolored water. (US only)	71406	4-32 oz
 Aqua Balance™ Maxibrom	Premium quality bromine is pH balanced. (EPA Registration #66397-1-1677) (US only)	17408	1-25 lb





Pool Parts Program

Your Ecolab Service Specialist can order pool equipment and accessories for your property and have it delivered on-site.

Test Kits

We can supply your facility with health-department approved test kits.


Aqua Balance™ Products (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Aqua Balance™ Poly 60	Controls algae growth. (EPA Registration #7124-77-1677) (US only)	17405	4-32 oz
 Aqua Balance™ Tinyl Foam	Removes grease, dirt and oils from vinyl pool walls and tile. (US only)	17407	4-32 oz





SeaKlear™ Products

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 SeaKlear™ Natural Clarifier	The all-natural way to remove dirt and oil from cloudy water. Enhances filter performance, safe to use with all sanitizer types and cannot be overdosed and over-mixed.	17572	4-32 oz
 SeaKlear™ 90-Day Algae Prevention & Remover for Pools	Kills and prevents all types of algae in water for 90 days, including blue-green, black and mustard. (EPA Registration #72083-1) (US only)	17573	4-32 oz
 SeaKlear™ Phosphate Remover	Removes up to 3,000 ppb of phosphates per 20,000 gallons of water without pool shutdown. (US only)	17571	4-32 oz
 SeaKlear™ Filter Cleaner and Degreaser	Restores filter efficiency and extends filter life by keeping filter elements clean and grease-free. (US only)	17561	4-32 oz
 SeaKlear™ PRS	Enhances filtration by trapping Cryptosporidium, E. coli, Giardia and algae with regular application. (US only)	17563	4-32 oz
 SeaKlear™ Self-Floccing Defoamer	Removes the source of foaming and scum without fouling your filter. (US only)	17562	4-32 oz



Filter Cartridges

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 ECL-T010	Three stage filtration technology utilizing hollow carbon, rated for 0.5 um particulate reduction, chlorine reduction of 50,000 gallons at 5.0 gpm and chloramine reduction of 8,000 gallons at 2.0 gpm.	9320-2198	10 in
 ECO-T014	Three stage filtration technology utilizing hollow carbon, rated for 0.5 um particulate reduction, chlorine reduction of 75,000 gallons at 7.0 gpm and chloramine reduction of 15,000 gallons at 3.7 gpm.	9320-2207	14 in
 ECO-T020	Three stage filtration technology utilizing hollow carbon, rated for 0.5 um particulate reduction, chlorine reduction of 100,000 gallons at 10.0 gpm and chloramine reduction of 22,500 gallons at 3.0 gpm.	9320-2250	20 in
 ECO-350	Three stage filtration technology utilizing hollow carbon, rated for chlorine reduction of 20,000 gallons at 1.7 gpm.	9320-2350	10 in
 ECO-315T05	Wet molded gradient density carbon block rated for 5 um particulate reduction and chlorine reduction of 6,000 gallons at 3.0 gpm. Starter Kit needed for install.	9320-1037 9320-1035 (installation kit)	16.1 in




Scale Control Filters

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 ECO-T010S	Three stage filtration technology utilizing hollow carbon, rated for 0.5 um particulate reduction, chlorine reduction of 50,000 gallons at 5.0 gpm and chloramine reduction of 8,000 gallons at 2.0 gpm.	9320-2256	10 in
 ECO-T014S5	Three stage filtration technology utilizing hollow carbon, rated for 0.5 um particulate reduction, chlorine reduction of 75,000 gallons at 7.0 gpm and chloramine reduction of 15,000 gallons at 3.7 gpm.	9320-2258	14 in
 ECO-T020S5	Three stage filtration technology utilizing hollow carbon, rated for 0.5 um particulate reduction, chlorine reduction of 100,000 gallons at 10.0 gpm and chloramine reduction of 22,500 gallons at 3.0 gpm.	9320-2259	20 in
 ECO-315T05S	Wet molded gradient density carbon block rated for 5 um particulate reduction and chlorine reduction of 6,000 gallons at 3.0 gpm. Starter Kit needed for install.	9320-1042 9320-1035 (Starter kit)	16.1 in 1 kit




Scale Control Filters (CONTINUED)

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 ECO-315CAT	Ion exchange (softening) capsule rated for 2,200 grains at 1.0 gpm. Starter Kit needed.	9320-1033 9320-1035 (Starter kit)	16.1 in 1 kit
 ECO SED10S CARTRIDGE	Sediment filter rated for 5 um particulate reduction at 5 gpm.	9320-2260	1 each


Equipment Systems

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 ECO-SINGLE	Single head filtration system.	9320-2266	1 each
 ECO-TWIN	Twin filtration system.	9321-2267	1 each
 ECO-SLF	Use with ECO-350 and ECO-350S.	9320-2270	1 each


Water Softeners

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 ECO-HTSS150	POU system, 8" x 17" tank size, 17" x 15" x 24" system size, 12" x 16" x 20" brine drum dimensions, 0.04 ft ³ resin volume per tank, 11 min/14 gal regeneration time/volume, 10.5-14.8 gpm flow range (15/30 psig).	9320-2237	1 each
 HWS-FLT	POU system, 8" x 17" tank size, 17" x 15" x 24" system size, 12" x 16" x 20" brine drum dimensions, 0.04 ft ³ , 11 min/14 gal regeneration time/volume, 10.5-14.8 gpm flow range (15/30 psig). (US only)	9320-2170	1 each
 EST-CP150	Whole house system, 10" x 54" tank size, 21" x 12" x 60" system size, 24" x 24" x 40" brine drum dimensions, 1.5 ft ³ resin volume per tank, 90 min/102 gal regeneration time/volume, 21-31 gpm flow range (15/30 psig). (US only)	9320-2148	1 each




Water Softeners


PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 EST-CP250	Whole house system, 13" x 54" tank size, 27" x 13" x 60" system size, 24" x 24" x 40" brine drum dimensions, 2.5 ft ³ resin volume per tank, 90 min/142 gal regeneration time/volume, 28-40 gpm flow range (15/30 psig). (US only)	9320-2149	1 each

Hot Water Filtration

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Hot Water Filtration Products	Stainless steel construction meant for laundry or warewashing, appropriate for temperatures under 275° F. Order cartridges separately.		
	Max flow rate 10, 250 psi dynamic pressure Replacement cartridge	9320-3248 9320-1023	1 each 1 cartridge
	Max flow rate 15, 250 psi dynamic pressure Replacement cartridge	9320-3255 9320-1024	1 each 1 cartridge
	Max flow rate 28, 125 psi dynamic pressure Replacement cartridge	9320-3263 9320-1023	1 each 1 cartridge
	Max flow rate 56, 125 psi dynamic pressure Replacement cartridge	9320-3271 9320-1024	1 each 1 cartridge

Optional Parts & Install Kits

PRODUCT	DESCRIPTION	PRODUCT CODE	PACK SIZE
 Braided Hose Install Kit	3/4" or 1/2" copper pipe.	9320-1094 9320-1095	1/2 in 3/8 in
 Bypass Plug	Bypass plug.	9320-2272	1 each
 Capsule Starter Kit	Mounting bracket for initial install of capsule filter.	9320-1035	1 kit

 **Replacement Cartridges**
Other legacy Ecolab and common competitive filter replacements are also available. Contact your Ecolab Representative for more information.

Unparalleled sales and service you can trust.



PREPARED

- Around-the-clock emergency service coverage
- Call 1-800-35-CLEAN 24/7/365

PROACTIVE

- Regularly scheduled service, sales and support
- On-site inspections for spot-free results
- Ongoing, on-site employee training

PROFESSIONAL

- Technical support from world-class chemists, microbiologists and engineers
- Optimize cleaning and disinfectant procedures
- Minimize energy and water costs