



## TSUNAMI™ 100: EXTENDING SHELF-LIFE

A Solution Story

Tsunami™ 100 helped a Southern **fresh-cut processor** extend shelf-life on fresh-cut greens.

### CHALLENGE

The processor was using acidified chlorine in their flume waters at 20 ppm active chlorine and pH of 6.7 and was seeking to **extend the shelf-life** of their product.

### SOLUTION

Replace acidified chlorine with **Tsunami 100** in the flume waters.

### RESULTS

REDUCED  
MICROBIAL  
POPULATIONS



SHELF-LIFE  
INCREASE  
with Tsunami 100




35%

for Collards  
& Kale



12%

for Turnip &  
Mustard Greens



To find out more, contact your Ecolab representative or call 1-800-392-3392.